



AUTOCLAVE



FOR SINGLE / MULTIPLE DRUMS



(ctrl+click on image)

Electrically operated.

Safety valve for release of steam pressure beyond 16 psi

Pressure gauge

Fully stainless steel body

Suitable for Drum Size: 9" x 9"

Also Available in Custom Size for 2 or 3 Drums

Specifications:

Construction	Fully Stainless Steel
Power Rating – Heater	1000 Watts to 2000 Watts
Pressure Control	By release of steam
Operation	Manual
Operating Voltage	230 Volts +/- 10%, 50 Hz, 1Phase
Pressure Indication	By Pressure Gauge.

DESIGN & SPECS SUBJECT TO CHANGE WITHOUT NOTICE

Manufactured By:

PUNEET INDUSTRIES

224, Sagar Ratan Industrial Premises, Plot D-265, TTC Industrial Area
Thane - Belapur Road, Turbhe MIDC, Navi Mumbai 400 705, INDIA.
+9198204 30094 puneetind@hotmail.com www.puneetind.com